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- Pawankar R., Canonica GW, Holgate S.T., Lockey R.F., WAO White Book on Allergy, 2011.
- Samolinski B., Raciborski F., Epidemiology of allergic diseases in Poland - ECAP study, *Allergy*, 2008;63:626-7.
- Penichet M.L., Jensen-Jarolim E., Cancer and IgE Introducing the Concept of Allergo-Oncology, Chapter 2: Humana Press, 2010.
- Leffell M., Donnenberg A., Rose N., Handbook of Human Immunology 1 ed.: Boca Raton FL: CRC Press.
- Kelso J.M., Sodhi N., Gosselin V.A., Yunginger J.W., Diagnostic performance characteristics of the standard Phadebas RAST, modified RAST, and Pharmacia CAP system versus skin testing, *Ann Allergy*, 1991 Nov;67(5):511-4.
- Pastorello E.A., Incorvaia C., Ortolani C., Bonini S., Canonica G.W., Romagnani S. et al., Studies on the relationship between the level of specific IgE antibodies and the clinical expression of allergy: I. Definition of levels distinguishing patients with symptomatic from patients with asymptomatic allergy to common aeroallergens, *J Allergy Clin Immunol*, 1995 Nov;96(5 Pt 1):580-7.
- Lavins B.J., Dolen W.K., Nelson H.S., Weber R.W., Use of standardized and conventional allergen extracts in prick skin testing, *J Allergy Clin Immunol*, 1992 Mar;89(3):658-66.
- Hahtela T., Jaakonmaki I., Relationship of allergen-specific IgE antibodies, skin prick tests and allergic disorders in unselected adolescents, *Allergy*, 1981 May;36(4):251-6.
- Chong Neto H.J., Rosario N.A., Studying specific IgE: in vivo or in vitro, *Allergol Immunopathol (Madr)*, 2009 Jan-Feb;37(1):31-5.
- Chapman M.D., Ferreira F., Villalba M., Cromwell O., Bryan D., Becker W.M. et al., The European Union CREATE project: a model for international standardization of allergy diagnostics and vaccines, *J Allergy Clin Immunol*, 2008 Nov;122(5):882-9 e2.
- Malandain H., IgE-reactive carbohydrate epitopes--classification, cross-reactivity, and clinical impact, *Eur Ann Allergy Clin Immunol*, 2005 Apr;37(4):122-8.
- Malandain H., IgE-reactive carbohydrate epitopes--classification, cross-reactivity, and clinical impact (2nd part), *Eur Ann Allergy Clin Immunol*, 2005 Sep;37(7):247-56.
- Foetisch K., Westphal S., Lauer I., Retzek M., Altmann F., Kolarich D. et al., Biological activity of IgE specific for cross-reactive carbohydrate determinants, *J Allergy Clin Immunol*, 2003 Apr;111(4):889-96.
- Vieths S., Luttkopf D., Reindl J., Anliker M.D., Wuthrich B., Ballmer-Weber B.K., Allergens in celery and zucchini, *Allergy*, 2002;57 Suppl 72:100-5.
- Sten E., Stahl Skov P., Andersen S.B., Torp A.M., Olesen A., Bindslev-Jensen U. et al., Allergenic components of a novel food, Micronesian nut Nangai (*Canarium indicum*), shows IgE cross-reactivity in pollen allergic patients, *Allergy*, 2002 May;57(5):398-404.
- Mari A., IgE to cross-reactive carbohydrate determinants: analysis of the distribution and appraisal of the in vivo and in vitro reactivity, *Int Arch Allergy Immunol*, 2002 Dec;129(4):286-95.
- van Ree R., Aalberse R.C., Demonstration of carbohydrate-specific immunoglobulin G4 antibodies in sera of patients receiving grass pollen immunotherapy, *Int Arch Allergy Immunol*, 1995 Feb;106(2):146-8.
- Iacovacci P., Afferni C., Butteroni C., Pironi L., Puggioni E.M., Orlandi A. et al., Comparison between the native glycosylated and the recombinant Cup a1 allergen: role of carbohydrates in the histamine release from basophils, *Clin Exp Allergy*, 2002 Nov;32(11):1620-7.
- Faye L., Gomord V., Fitchette-Laine A.C., Chrispeels M.J., Affinity purification of antibodies specific for Asn-linked glycans containing alpha 1-->3 fucose or beta 1-->2 xylose, *Anal Biochem*, 1993 Feb 15;209(1):104-8.
- Lauriere M., Lauriere C., Chrispeels M.J., Johnson K.D., Sturm A., Characterization of a xylose-specific antiserum that reacts with the complex asparagine-linked glycans of extracellular and vacuolar glycoproteins, *Plant Physiol*, 1989 Jul;90(3):1182-8.
- Ramirez-Soto D., Poretz R.D., The (1-->3)-linked alpha-L-fucosyl group of the N-glycans of the *Wistaria floribunda* lectins is recognized by a rabbit anti-serum, *Carbohydr Res*, 1991 Jun 25;213:27-36.
- van der Veen M.J., van Ree R., Aalberse R.C., Akkerdaas J., Koppelman S.J., Jansen H.M. et al., Poor biologic activity of cross-reactive IgE directed to carbohydrate determinants of glycoproteins, *J Allergy Clin Immunol*, 1997 Sep;100(3):327-34.
- Altmann F., The role of protein glycosylation in allergy, *Int Arch Allergy Immunol*, 2007;142(2):99-115.